

Mixing Directions - Sterilex Liquid Ultra

DIRECTIONS FOR USE:

1. Initial Start-Up treatment: treat the DUWL on three consecutive nights.
2. Routine treatment: treat DUWL one night/week after initial start-up treatment.
3. If two or more weekly treatments are skipped, repeat initial start-up treatment.

HOW TO APPLY STERILEX LIQUID ULTRA:

1. Add one bottle (3 ounces) of Solution I and one bottle (3 ounces) of Solution II to an empty Sterilex measuring cup or other empty container and stir. (Note: Solution I and II must be mixed together prior to treatment, and the mixture must be used within 24 hours.) Pour the appropriate amount of combined mixture into an empty external dental unit water bottle. Do not use PETE bottles (see note below).

2. Purge lines. Run Sterilex Liquid Ultra mixture through the system until the pink solution appears at the end of the A/W syringes and handpiece lines. **Always remove handpiece. It may be advisable to remove coupler depending upon the type of coupler. Contact manufacturer for specific recommendations.**

3. Allow the Sterilex Liquid Ultra mixture to sit in the unit overnight (place the ends of waterlines into a sink or cuspidor in case any pink solution drips overnight).

4. At the beginning of the next workday, purge lines of Sterilex Liquid Ultra mixture. Rinse the external water bottle with water.

5. Fill external bottle with water. Flush each line (A/W syringe, handpiece lines) until the pink solution is no longer visible.

SAFETY INFORMATION:

Read instructions thoroughly. Avoid contact with skin and eyes. If product gets in eyes or on skin flush thoroughly. Do not ingest. If swallowed, drink large quantities of water. Get immediate medical attention.

Keep out of children's reach. Store at room temperature out of direct sunlight, do not store above 30°C (or 86°F).

Note: Do not use transparent polyethylene terephthalate (PETE) bottles. We recommend either a high density polyethylene bottle with a minimum thickness of 0.08 inches or a high/low density polyethylene blend with a minimum thickness of 0.14 inches.

Control of bacterial contamination in the DUWL is affected by the quality of incoming water. Sterile water is recommended by OSAP.