

Mixing Directions - Powder

DIRECTIONS FOR USE:

1. Initial start-up treatment: treat the DUWL on three consecutive nights.
2. Routine treatment: treat DUWL one night/week after initial start-up treatment.
3. If two or more weekly treatments are skipped, repeat initial start-up treatment.

HOW TO APPLY STERILEX ULTRA:

1. At the end of the workday, add eight (8) ounces of **HOT** water to a Sterilex measuring cup or other empty container.
2. Add one packet of Sterilex Ultra to the **HOT** water. Use stirrer to mix or swirl to dissolve. Be sure all powder is dissolved. Pour appropriate amount of pink liquid into an empty external dental unit water bottle. Solution must be used within 24 hours of mixing. Do not use PETE bottles (see below note for bottle requirements).
3. Purge lines. Run Sterilex Ultra solution through the system until the pink solution appears at the end of the a/w syringe and handpiece lines. Always remove handpiece. **It may be advisable to remove the coupler depending upon the type of coupler. Contact manufacturer for specific recommendations.**
4. Allow the Sterilex Ultra solution to sit in the unit overnight (place the ends of waterlines into a sink or cuspidor in case any pink solution drips overnight).
5. At the beginning of the next workday, purge lines of Sterilex Ultra solution. Rinse the external water bottle with hot water.
6. Fill bottle with hot water. Flush each line (a/w syringe, handpiece lines) until the bottle is empty.

SAFETY INFORMATION:

Read instructions thoroughly. Avoid contact with skin and eyes. If product gets in eyes or on skin flush thoroughly. Do not ingest. If swallowed, drink large quantities of water. Get immediate medical attention.

Keep out of children's reach. Store at room temperature out of direct sunlight, do not store above 30°C (or 86°F).

Note: Do not use transparent polyethylene terephthalate (PETE) bottles. We recommend either a high density polyethylene bottle with a minimum thickness of 0.08 inches or a high/low density polyethylene blend with a minimum thickness of 0.14 inches.

Control of bacterial contamination in the DUWL is affected by the quality of incoming water. Sterile water is recommended by OSAP.